



• xFuse Chips DBM Powder:

(a mix of demineralized bone and cancellous chips) available in a particle sizes ranging from 125-850um and in different volumes formulations. xFuse DBM Powder provide both an osteoinductive component as well as an osteoconductive potential.

SIZE 125-850um	PRODUCT CODE
lcc	ALLO-DBM-01
200	ALLO-DBM-02
5cc	ALLO-DBM-05
10cc	ALLO-DBM-10

Product Features

xFuse Chips are stored in ambient temperature featuring "off the shelf technology"

xFuse Chips have a Sterility Assurance Level of 10⁻⁶ (SAL) Provide an onsteoinductive / osteoconductive matrix potential Easy to mix with autograft or xFuse Fiber Boats and Putties as a graft extender.

100% Allograft.

info@allotekinc.com
4040 Broadway St.
Suite 240-308
San Antonio, Texas, 78209



CALL CUSTOMER SERVICE



Allotekinc. **XFuse Cancellous Powder**



• xFuse Chips Cancellous Powder:

are available in particle size ranges (.5mm - 1mm) and four volumes and are a mix of crushed cancellous and crushed cortical chips ideal for dental and maxillofacial needs.

SIZE .5 - 1mm	PRODUCT CODE
.5cc	ALLO-CB-00
100	ALLO-CB-01
200	ALLO-CB-02
5cc	ALLO-CB-05



Product Features

xFuse Chips are stored in ambient temperature featuring "off the shelf technology"

xFuse Chips have a Sterility Assurance Level of 10⁻⁶ (SAL) Provide an onsteoinductive / osteoconductive matrix potential Easy to mix with autograft or xFuse Fiber Boats and Putties as a graft extender.

100% Allograft.

info@allotekinc.com 4040 Broadway St. Suite 240-308 San Antonio, Texas, 78209



CALL CUSTOMER SERVICE

1 877 2556835 1 877ALLOTEK





xFuse Chips Cortical Ground Filler

<1mm: are a mix of crushed cancellous and crushed cortical chips ideal for dental and maxillofacial needs or for packing voids.

SIZE <1mm	PRODUCT CODE
.5cc	ALLO-GB-00G
lcc	ALLO-GB-01G
2.5cc	ALLO-GB-02G
5cc	ALLO-GB-05G
10cc	ALLO-GB-10G
20cc	ALLO-GB-20G
30cc	ALLO-GB-30G

Product Features

xFuse Chips are stored in ambient temperature featuring "off the shelf technology"

xFuse Chips have a Sterility Assurance Level of 10⁻⁶ (SAL) Provide an onsteoinductive / osteoconductive matrix potential Easy to mix with autograft or xFuse Fiber Boats and Putties as a graft extender.

100% Allograft.

info@allotekinc.com
4040 Broadway St.
Suite 240-308
San Antonio, Texas, 78209



CALL CUSTOMER SERVICE



Allotek inc. **Fuse Chips Cancellous



• xFuse Chips Cancellous:

are available in two particle size ranges (1-4mm and 4-10mm) and five standard volumes and are a standard in a wide range of surgical procedures.

SIZE 1-4mm	PRODUCT CODE
5cc	ALLO-CB-05G
15cc	ALLO-CB-15G
30cc	ALLO-CB-30G
60cc	ALLO-CB-60G

SIZE 4-10mm	PRODUCT CODE
5cc	ALLO-CB410-05G
10cc	ALLO-CB410-10G
15cc	ALLO-CB410-15G
30cc	ALLO-CB410-30G
60cc	ALLO-CB410-60G



Size 1-4mm



Size 4-10mm

info@allotekinc.com 4040 Broadway St. Suite 240-308 San Antonio, Texas, 78209



CALL CUSTOMER SERVICE

